Reg. Numero / Reg. Number	MDR 00054-A	Revisione / Revision	0
Primo rilascio / First issue date	2024-05-31	Valido da / Valid from	2024-05-31
Scadenza / Valid until	2029-05-30	Ultima modifica / Last change date	2024-05-31
Data Scadenza Precedente / Previous Expiry			

Pagina 1 di 6
Page 1 of 6

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Kiwa Cermet Italia certifies that, on the basis of the assessment carried out, the Quality Management System of the Organization:

MEDITALIA INDUSTRIALE S.r.l.

Operatore economico / Economic operator: Fabbricante / Manufacturer

SRN: IT-MF-000032344

Sede Legale e Operativa / Legal and Operational Headquarters

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Is in compliance with the applicable requirements of Regulation (EU) 2017/745, Annex IX Chapter I and III, for the following devices types:

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

Solutran[®]hemo BioNest Transfer Bag configuration

Field of use

System for the production and storage of blood components from umbilical cord blood.

Warnings

1. Bring the temperature of the bags between 15 ° C and 30 ° C before use.
2. Heating in a microwave oven is not recommended.
3. Follow aseptic procedures.
4. Avoid contact with sharp objects.
5. Inspect the packaging and the bags and DO NOT use if there are any defects or anomalies.
6. DO NOT use if the primary packaging is open or damaged.
7. Check the expiration date before use.
8. DO NOT vent.
9. DO NOT add drugs.
10. Single use device. DO NOT reuse. Risks related to re-use: infection/sepsis, coagulation, loss of efficacy.
11. DO NOT re-sterilize. Re-sterilization could cause cross-contamination by patient or operator

12. Follow validated procedures for the preparation of blood components. In case of blood fractionation by automatic separator, follow the user manual of the separator.
13. Before use, blood components must be evaluated according to the parameters required by applicable laws or standards.
14. This medical device is intended for use by qualified medical personnel only.

Instructions for Use

1. Follow aseptic procedures.
2. Seal tubes close to caps no. 4 from BC bags.
3. Clamp tubes connecting bags A to B and B to C.
4. Aseptically connect whole CB unit to BioNest bag A by sterile connection device.
5. Transfer unit identification data on BioNest ABC bags according to local procedures.
6. Transfer whole CB from original CB collection bag into BioNest bag A. Separate original CB container.
7. Use flexible protective envelope and solid cylinder adaptor (provided separately) suitable for the centrifuge buckets. Insert ABC bags into flexible protective envelope and solid cylinder adaptor.
8. Centrifuge ABC bags at low speed (e.g. 200-240 g for 10 minutes) to separate the CBPRP (top of bag A), from the CBRBC (bottom of bag A).
9. Carefully remove ABC bags from solid cylinder adaptor and flexible protective envelope, set bag A into manual plasma extractor and slowly transfer CBPRP into bag B.
10. Repeat step 7.
11. Centrifuge ABC bags at high speed (e.g. 2000 g for 15 minutes) to separate the low-hematocrit CBRBC and the CBPRP obtained in step 8 into high-hematocrit CBRBC (bottom of bag A) + platelet poor plasma (CBPPP, top of bag A) and platelet concentrate (CBPC, bottom of bag B) + CBPPP (top of bag B).
12. Carefully remove ABC bags from solid cylinder adaptor and flexible protective envelope, set bag A into manual plasma extractor and slowly transfer CBPPP from top of bag A into bag C. Pay attention not to disturb the bottom layer of platelets concentrated in bag B.
13. Separate bag A from bags BC, store CBRBC at 2-6 °C and gently set bag B into manual plasma extractor.
14. Transfer supernatant CBPPP from bag B into bag C except the volume necessary to achieve a pre-defined platelet concentration in the CBPC (e.g. 800-1200 x 10⁹/L).
15. Separate bags B and C, resuspend the platelets in the CBPC. Store CBPC and CBPPP according to local procedures. Bags B and C can be cryopreserved at temperatures between -25 °C and -80 °C.
16. For the preparation of cord blood platelet gel (CBPG), remove the open lock cap from cap n 10 and, using sterile alcohol wipe, swab the surface of valve and let air dry; connect cap 10 of bag B to cap 3 of bag D "Solutran@hemo BioNest Transfer Bag" and transfer the platelet concentrate into bag D, which can be frozen at temperatures between -25 and -80 °C. Seal connecting tube and discard bag B.
17. Prepare platelet gel by aseptically injecting through cap no. 2 a volume of a 10% calcium gluconate solution corresponding to one third of the volume of the platelet concentrate into bag D (after defrosting). The gel is formed at room temperature in about 30-40 minutes. The extraction of the gel and its topical application is facilitated by the tear opening along the lateral engraving of bag D.

Codes

4911002200 Solutran@hemo BioNest System for UCB

4211000000 Solutran[®]hemo BioNest Transfer Bag

Storage conditions

Store in a cool, dry place, away from heat and direct sunlight.

Instructions for disposal:

Follow norms and regulations for potentially contaminated biological products

Symbol	Meaning	Symbol	Meaning
	Avoid direct exposure to sun light and heat sources		Manufacturer
	Product Code		Lot
	Expiration date		Manufacturing date
	Sterilized by ethylene oxide		Single use
	Do not use if package is damaged or deteriorated		Do not vent
	Pyrogen free fluid path		Do not re-sterilize
	Phthalate (DEHP) free		Contains phthalates (DEHP)
	Instructions for use		Whole blood bag
	Plasma bag		Platelet bag
	Do not use if bags show evidence of damage		Transfer bag
	Compliant with Regulation (EU) 2017/745.		

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

SAGM sterile bag Configuration

Field of use

System for the production and storage of leukoreduced cord blood red blood cells (CBRBC).

Warnings

1. Bring the temperature of the BioNest EF system between 15 °C and 30 °C before use.
2. Heating in a microwave oven is not recommended.
3. Follow aseptic procedures.
4. Avoid contact with sharp objects.
5. Inspect the packaging and the bags and DO NOT use if there are any defects or anomalies.
6. DO NOT use if the primary packaging is open or damaged.
7. Check the expiration date before use.
8. DO NOT vent.
9. DO NOT add drugs.
10. Single use device. DO NOT reuse. Risks related to re-use: infection/sepsis, coagulation, loss of efficacy.
11. DO NOT re-sterilize. Re-sterilization could cause cross-contamination by patient or operator
12. Follow validated procedures for the preparation of blood components. In case of blood fractionation by automatic separator, follow the user manual of the separator.
13. Before use, blood components must be evaluated according to the parameters required by applicable laws or standards.
14. This medical device is intended for use by qualified medical personnel only.

Instructions for Use

1. Follow aseptic procedures.
2. Remove BioNest EF from the container.
3. Clamp all tube segments.
4. Aseptically connect high-hematocrit (e.g. 80-90%) cord blood red blood cells (CBRBC) to BioNest EF system.
5. Aseptically connect a bag containing sterile SAGM bag to segment no. 12. The Meditalia Industriale product code 710000S000 can be used for this operation; in this case, remove the SAGM bag from over pouch, remove immediately open lock cap from cap n° 4 of SAGM bag and connect immediately aseptically to segment n° 12 of BioNest EF system.
6. Transfer pre-defined SAGM volume to high-hematocrit CBRBC (e.g. 0.5 vol SAGM to 1 vol CBRBC) to achieve desired CBRBC hematocrit (e.g. 55-60%). Clamp segment no. 12 and open clamp from segment no. 13.
7. Filter by gravity.

8. Collect leukoreduced CBRBC in storage bag.
9. Store leukoreduced CBRBC at 2-6 °C according to local procedures.

Codes

4310202000 Solutran[®] hemo BioNest EF

710000S000 SAGM sterile bag

Storage conditions

Store in a cool, dry place, away from heat and direct sunlight.

Instructions for disposal:

Follow norms and regulations for potentially contaminated biological products

Symbol	Meaning	Symbol	Meaning
	Avoid direct exposure to sun light and heat sources		Manufacturer
	Product Code		Lot
	Expiration date		Manufacturing date
	Sterilized by ethylene oxyde		Steam sterilization
	Do not use if package is damaged or deteriorated		Single use
	Pyrogen free fluid path		Do not vent
	Do not re-sterilize		Phthalate (DEHP) free
	Contains phthalates (DEHP)		Instructions for use
	Do not use if bags show evidence of damage		Transfer bag (collection bag of leukoreduced cbrbc)
	Compliant with Regulation (EU) 2017/745.		

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
Field of use

Bag system for the preparation of multiple aliquots of platelet gel from adult donors.

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Warnings

1. Let the bag temperature reach 15 – 30 °C before use.
2. Heating in a micro-oven device is not recommended.
3. Follow aseptic procedures.
4. Avoid contact with sharp objects.
5. Inspect the packaging and the bags and do not use them if defects and anomalies are detected.
6. Do not use if the primary packaging is open or damaged.
7. Check expiration date before use.
8. Do not vent.
9. Single use device. Do not re-use. Re-use risks: infection/sepsis, coagulation, loss of efficacy.
10. Do not re-sterilize. Re-sterilization could cause patient or operator cross-contamination.
11. Follow validated procedures for the preparation of blood components. In case of blood fractionation by automatic separator, follow the user manual of the separator.
12. Before use, blood components must be evaluated according to the parameters required by applicable laws or standards. (In Italy, Decree. 2/11/15).
13. This medical device is intended for use by qualified medical personnel only.
14. Be sure that product handling during its use is produced in adequate or controlled environments (ex. Environments under periodic controls or sterile flow cabinets) and that any connections is performed with sterile connectors or in sterile flow cabinets and by qualified personnel.

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Solutran[®]hemo BioNest System for adult platelet donors

Instructions for use

1. Follow aseptic procedures.
2. Connect the system to the bag containing the blood component to be divided into aliquots.
3. Transfer aliquots as appropriate into the transfer bags.
4. Seal the connecting tubes.
5. Complete the labelling procedure.
6. Separate the bags and discard the original blood component empty bag.

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Instructions for disposal:

Follow norms and regulations for potentially contaminated biological products.

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Codici applicabili

4511002400 Quadruple bag system for platelet storage 150 mL
4511007400 Quadruple bag system for platelet storage 400 mL
4511003400 Quadruple bag system for plasma storage 200 mL
4511007401 Quadruple bag system for plasma storage 400 mL
4511003401 Quadruple bag system for red blood cell storage 200 mL
4511007402 Quadruple bag system for red blood cell storage 400 mL
4911003600 Sextuple bag system for red blood cell storage 200 mL
4911007600 Sextuple bag system for red blood cell storage 400 mL

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Codici applicabili

4111003000 Single bag for blood components storage 200 mL
4111007000 Single bag for blood components storage 400 mL
4111011000 Single bag for blood components storage 600 mL
4111011001 Single bag for red blood cell storage 600 mL
4111019100 Single DEHP free bag for platelet storage, one connection tube, 1000 mL
4211020100 Single DEHP free bag for platelet storage, two connection tubes, 1300 mL
4110020000 Single DEHP free bag for platelet storage, one connection tube, 1300 mL
4311007100 Single DEHP free bag for platelet storage with sampling bag, one connection tube, 400 mL
4311011100 Single DEHP free bag for platelet storage with sampling bag, one connection tube, 600 mL
4311019100 Single DEHP free bag for platelet storage with sampling bag, one connection tube, 1000 mL
4311020100 Single DEHP free bag for platelet storage with sampling bag, one connection tube, 1300 mL

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